

Qualitative Fit Test Kits

The biggest contributor to reduced respiratory protection is poor fit. Conducting a fit test to ensure that a respirator with a tight fitting facepiece provides an adequate seal to the wearer's face has long been considered best practice as part of a general respiratory protection program.

Reflecting this, a fit test for wearers of tight fitting respirators is now required for compliance to AS/NZS1715:2009.

The 3M Qualitative Fit Test apparatus FT-10 and FT-30 are designed for fit testing disposable and half face reusable respirators fitted with particulate or combination gas/vapour and particulate filters.

Features and Benefits

- Offers a fast and easy method for performing qualitative fit testing
- No calibration of equipment required
- No modification of facepiece required
- Test operator can be self-taught/instructed
- Each kit contains a hood and collar assembly, two nebulisers, sensitivity solution, test solution and detailed instructions
- Replacement solutions and accessories can be purchased separately
- Available with Sweet (FT-10) or Bitter (FT-30) solutions



For more information, please visit
www.limkimhai.com.sg/3m-sensitivity-fit-test-solutions/

Authorised Distributor:



Lim Kim Hai Electric Co (S) Pte Ltd
Lim Kim Hai Building,
53 Kallang Place, Singapore 339177
T (+65) 6292 3711 • (+65) 6490 5000
F (+65) 6297 0078
E customerservice@limkimhai.com.sg
W www.limkimhai.com.sg

Accessories

FT-13 Nebuliser
Packaging: 3/case

FT-14 Test Hood
Packaging: 2/pack;
5packs/case

3M™ Qualitative Fit Test Apparatus

FT-10 Qualitative Fit Test Apparatus Kit
- Sweet (Saccharin)

FT-30 Qualitative Fit Test Apparatus
Kit - Bitter (Bitrex)

Packaging FT-10 & FT-30: 1 kit/case

Sensitivity & Fit Test Solutions

FT-11 Sensitivity Solution (sweet)

FT-12 Fit Test Solution (sweet)

FT-31 Sensitivity Solution (bitter)

FT-32 Fit Test Solution (bitter)

*Packaging for all solution:
55ml bottle, 6 bottles/case*