



3M™ Quat Disinfectant Cleaner *(Ready to use)*

Product Details

- EPA registered disinfectant cleaner for use in hospitals
- Effective against a broad spectrum of bacteria, is virucidal, including HIV-1, the AIDS Virus and Hepatitis B Virus and Hepatitis C Virus, fungicidal and inhibits the growth of mold and mildew and their odors
- Proven “one-step” disinfectant / cleaner effective against a broad spectrum of bacteria, virus & fungi
- This product is a disinfectant that cleans, disinfects, deodorizes with a pleasant fragrance and rinse-free in one labour-saving step
- Applications: Washable, hard, nonporous surfaces of floors, walls, countertops, sinks, tables, chairs, garbage cans, telephones, bathtubs, urinals, toilet bowls, and other hard, nonporous surfaces made of metal, glazed ceramic, plastic and stainless steel
- Volumes Available: 650ml (Most Popular), 4 Litre, 20 Litre



Test Result

	Before Cleaning	After Cleaning
Floor	1397 RLU	103 RLU ↓
Door	339 RLU	140 RLU ↓

* Relative Light Unit (RLU) is a unit for measuring cleanliness by measuring the levels of Adenosine Triphosphate.

- ❖ The door is applied with Quat Disinfectant Cleaner using mist spray or blower machine and then air-dried (Treated surfaces can be wipe dry with a cloth, sponge, mop or allow to air dry)
- ❖ A Walk - Behind Scrubber and Dryer machine is used to disperse Quat Disinfectant Cleaner and then scrub and dry the floor.



Lim Kim Hai Electric Co (S) Pte Ltd
Lim Kim Hai Building,
53 Kallang Place, Singapore 339177
www.limkimhai.com.sg

Contact us today at solutions@limkimhai.com.sg

Hotline: +65 6490 5000